



Group Summary Report

ULTIMATE-D BILIRUBIN XXII

Lot To Date through Oct 2009

CONSTITUENT UNIT OF MEASURE METHOD PRINCIPLE INSTRUMENT/SYSTEM	M804561/L1						M804562/L2						M804563/L3					
	AVG. MEAN	SD OF MEAN	AVG S.D.	AVG. C.V.	NO. FLS	NO. LABS	AVG. MEAN	SD OF MEAN	AVG S.D.	AVG. C.V.	NO. FLS	NO. LABS	AVG. MEAN	SD OF MEAN	AVG. S.D.	AVG. C.V.	NO. FLS	NO. LABS
Bilirubin, Direct																		
mg/dL																		
Diazo-Other (J-G) w/B																		
BCI CX 3/4/5/7/9	0.729	0.028	0.054	7.4	38	38	3.611	0.057	0.077	2.1	18	18	5.304	0.080	0.125	2.4	39	39
BCI LX Systems	0.779	0.030	0.055	7.1	25	22	3.617	0.034	0.093	2.6	13	10	5.303	0.067	0.115	2.2	28	23
BCI DxC	0.773	0.037	0.055	7.1	110	83	3.599	0.055	0.099	2.8	54	42	5.291	0.063	0.119	2.3	106	84
All Instruments	0.764	0.039	0.055	7.1	173	143	3.604	0.053	0.094	2.6	85	70	5.296	0.068	0.120	2.3	173	146
All Method Principles	0.764	0.039	0.055	7.1	173	143	3.604	0.053	0.094	2.6	85	70	5.296	0.068	0.120	2.3	173	146
Bilirubin, Total																		
mg/dL																		
Diazo-Other(J-G) wo/B																		
BCI CX 3/4/5/7/9	1.38	0.11	0.10	7.3	13	13	5.45	0.06	0.14	2.5	8	8	7.99	0.07	0.23	2.9	12	12
BCI LX Systems	1.33	0.06	0.11	8.2	9	8	5.33		0.17	3.2	6	5	7.97	0.08	0.21	2.6	10	9
BCI DxC	1.31	0.07	0.11	8.3	46	35	5.33	0.08	0.16	3.0	25	19	7.98	0.11	0.21	2.7	44	35
All Instruments	1.32	0.08	0.11	8.1	68	56	5.36	0.09	0.16	2.9	39	32	7.98	0.10	0.22	2.7	66	56
Enz Bili Ox w/B																		
BCI DxC																		
All Instruments																		
All Method Principles	1.32	0.08	0.11	8.1	68	56	5.36	0.09	0.16	2.9	39	32	7.98	0.10	0.22	2.7	66	56



Group Summary Report

ULTIMATE-D BILIRUBIN XXII

Lot To Date through Oct 2009

CONSTITUENT UNIT OF MEASURE METHOD PRINCIPLE INSTRUMENT/SYSTEM	M804564/L4					
	AVG. MEAN	SD OF MEAN	AVG S.D.	AVG. C.V.	NO. FLS	NO. LABS
Bilirubin, Direct						
mg/dL						
Diazo-Other (J-G) w/B						
BCI CX 3/4/5/7/9	10.105		0.208	2.1	2	2
BCI LX Systems	10.206		0.128	1.3	2	2
BCI DxC	10.235		0.147	1.4	2	2
All Instruments	10.182		0.161	1.6	6	6
All Method Principles	10.182		0.161	1.6	6	6
Bilirubin, Total						
mg/dL						
Diazo-Other(J-G) wo/B						
BCI CX 3/4/5/7/9	16.70		0.27	1.6	3	3
BCI LX Systems	16.77	0.17	0.43	2.5	12	8
BCI DxC	16.71	0.24	0.36	2.2	23	19
All Instruments	16.73	0.21	0.37	2.2	38	30
Enz Bili Ox w/B						
BCI DxC	17.70		0.12	0.7	1	1
All Instruments	17.70		0.12	0.7	1	1
All Method Principles	16.75	0.26	0.37	2.2	39	31