



## Group Summary Report

ULTIMATE-D BILIRUBIN XXIII

Lot To Date through Oct 2009

CONSTITUENT UNIT OF MEASURE METHOD PRINCIPLE INSTRUMENT/SYSTEM	M812191/L1						M812192/L2						M812193/L3					
	AVG. MEAN	SD OF MEAN	AVG S.D.	AVG. C.V.	NO. FLS	NO. LABS	AVG. MEAN	SD OF MEAN	AVG S.D.	AVG. C.V.	NO. FLS	NO. LABS	AVG. MEAN	SD OF MEAN	AVG. S.D.	AVG. C.V.	NO. FLS	NO. LABS
<b>Bilirubin, Direct</b>																		
mg/dL																		
Diazo-Other (J-G) w/B																		
BCI CX 3/4/5/7/9	0.920		0.040	4.3	1	1	3.937		0.062	1.6	2	2	5.694		0.076	1.3	6	6
BCI LX Systems	1.000		0.034	3.4	1	1	3.954		0.089	2.3	1	1	5.686		0.077	1.4	3	2
BCI DxC	0.987	0.035	0.048	4.9	14	9	4.047	0.160	0.113	2.8	8	6	5.685	0.091	0.126	2.2	41	32
All Instruments	0.983	0.037	0.047	4.8	16	11	4.019	0.143	0.101	2.5	11	9	5.686	0.084	0.117	2.1	50	40
All Method Principles	0.983	0.037	0.047	4.8	16	11	4.019	0.143	0.101	2.5	11	9	5.686	0.084	0.117	2.1	50	40
<b>Bilirubin, Total</b>																		
mg/dL																		
Diazo-Other(J-G) wo/B																		
BCI CX 3/4/5/7/9							5.91		0.13	2.2	2	2	8.62		0.17	1.9	2	2
BCI LX Systems	1.60		0.15	9.2	1	1	6.00		0.14	2.4	1	1	8.73		0.27	3.1	2	1
BCI DxC	1.58	0.05	0.09	5.5	9	5	5.96		0.18	3.1	4	2	8.59	0.17	0.22	2.6	21	16
All Instruments	1.58	0.05	0.09	5.8	10	6	5.95	0.11	0.16	2.7	7	5	8.60	0.16	0.22	2.6	25	19
All Method Principles	1.58	0.05	0.09	5.8	10	6	5.95	0.11	0.16	2.7	7	5	8.60	0.16	0.22	2.6	25	19



**Group Summary Report**  
 ULTIMATE-D BILIRUBIN XXIII

Lot To Date through Oct 2009

CONSTITUENT UNIT OF MEASURE METHOD PRINCIPLE INSTRUMENT/SYSTEM	M812194/L4							
	AVG. MEAN	SD OF MEAN	AVG S.D.	AVG. C.V.	NO. FLS	NO. LABS		
Diazo-Other (J-G) w/B BCI CX 3/4/5/7/9 BCI LX Systems BCI DxC All Instruments  All Method Principles Diazo-Other(J-G) wo/B BCI CX 3/4/5/7/9 BCI LX Systems BCI DxC All Instruments  All Method Principles								